For microscopy

For macro photography

Ultra
33x

Micro
20x

Macro
10x

Macro Plus
5x

Video tutorial:

Experience our app:

Transform your smartphone into
a digital microscope!

blips
magnify the world

A product by

SMO
SMART MICRO OPTICS

SmartMicroOptics S.r.l.
Via G. di Cornigliano 6r
16152 Genova - Italy

new
LABKIT₂

www.blips.pro

    @blipslens

INSTRUCTIONS – new Blips LabKit 2

1. **Lenses** – Macro Plus (yellow), Macro (green), Micro (blue) and Ultra (red), provide the magnification of 5x, 10x, 20x and 33x, respectively. The Macro Plus and Macro lenses do not need the supporting phone stage and they can be used anytime for taking macro pictures. Micro and Ultra lenses need the supporting stage to provide their best performance, due to their high magnification rate.

2. **Lens on smartphone** – Stick the lens, well centered onto the phone camera (on the main camera, for multi-cam phones).

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 5. **Phone** – Lay a prepared slide on the stage over the light. Lay the phone on the U-shaped stands, with the camera aligned with sample and light source. Tilt the smartphone slightly, moving the lens away and closer to the sample, to find out if you have to move the specimen closer, or farther, to get it well in focus. Rotate the feet to raise or lower the sample accordingly. When you are close to the ideal focus rotate only one foot, for a fine adjustment. Check the alignment of the LED light.

6. **Your own sample** – Lay the slide with your own preparation on the hole of the stage. Important: keep the sample very thin. The transparent foils can be used as coverslips for your preparations. Cut them, if necessary.

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
 3. **Supporting stage preparation** - Stick the adhesive rubber on the U-shaped stands. Screw the rubber feet in the metallic plate for about 5-10 mm.

4. **Lighting** – Fix the plastic guide to the table with its adhesive pads (optional). Put the rubber feet of the slide holder inside the external holes of the plastic guide. Insert the light source in the central hole, with the led by the side of the small hole in the slide support. You will get a perfect alignment of the light source with the sample.

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 7. **Digitally zoom** the sample for a fine adjustment of focus and light alignment.

We suggest to lock the automatic focusing.

8. Find examples, information, videos and our free SMO app at www.blips.pro

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WARNING:
CHOKING HAZARD - Small parts.
Not for children under 3 years.