


Use DIPLE Lux like a Powerful magnifying lens for smartphone:




Turn-on the light:

When the external light is not enough to properly illuminate your sample, you can turn on two little LEDs under the DIPLE tile. With the switch ON/OFF you open and close the circuit.

With the + and - tactile switches you can increase or decrease the light intensity.


The light is off when you switch the to ON state. This is to avoid the batteries discharging in case of accidental switch. PRESS THE + (plus) BUTTON LONGER THAN 5-6 SECONDS TO INCREASE the intensity.

PRESS THE - (minus) BUTTON TO DECREASE the intensity.



Place the lens tile on the subject to observe. Keep DIPLE Lux in one hand, like a magnifying lens, and the phone on the other hand. With the camera-app on, move the phone toward the centre of the black rubber ring and get in contact with it. Slightly tilt the DIPLE Lux, or use the elevation screw in the tile to find the right working distance. Change the light intensity with the buttons + and -.

Use DIPLE Lux like a Powerful magnifying lens for smartphone:




Turn-on the light:

When the external light is not enough to properly illuminate your sample, you can turn on two little LEDs under the DIPLE tile. With the switch ON/OFF you open and close the circuit.

With the + and - tactile switches you can increase or decrease the light intensity.


The light is off when you switch the to ON state. This is to avoid the batteries discharging in case of accidental switch. PRESS THE + (plus) BUTTON LONGER THAN 5-6 SECONDS TO INCREASE the intensity.

PRESS THE - (minus) BUTTON TO DECREASE the intensity.



Place the lens tile on the subject to observe. Keep DIPLE Lux in one hand, like a magnifying lens, and the phone on the other hand. With the camera-app on, move the phone toward the centre of the black rubber ring and get in contact with it. Slightly tilt the DIPLE Lux, or use the elevation screw in the tile to find the right working distance. Change the light intensity with the buttons + and -.

xn7




xn7



DIPLE Lux

The powerful magnifying lens for your phone

Thanks for purchasing DIPLE®, a product designed and manufactured by SmartMicroOptics Sri. With your smartphone and DIPLE® LUX you can have a very powerful digital magnifying lens, anytime. You do not have to mount or stick anything on the phone; just move the objective lens of your phone toward the lens of the DIPLE tile and get in contact with it. You can switch on the light of DIPLE Lux toward your sample, if needed. You can tune the lens-subject working distance tilting the tile, like a classical magnifying lens, or using the provided elevation screw.




www.smartmicrooptics.com

DIPLE Lux

The powerful magnifying lens for your phone

Thanks for purchasing DIPLE®, a product designed and manufactured by SmartMicroOptics Sri. With your smartphone and DIPLE® LUX you can have a very powerful digital magnifying lens, anytime. You do not have to mount or stick anything on the phone; just move the objective lens of your phone toward the lens of the DIPLE tile and get in contact with it. You can switch on the light of DIPLE Lux toward your sample, if needed. You can tune the lens-subject working distance tilting the tile, like a classical magnifying lens, or using the provided elevation screw.




www.smartmicrooptics.com



DIPLE Lux:

The Diple Lux objective lens is connected to the handle. Inside the handle there is the circuit for the lighting of the sample



Batteries:




Two CR2032 (1.5V) are inside.

For changing/mounting the batteries use a screwdriver to release the three screws on the plastic shell.

Technical Specifications*:

Working distance (WD): 2.8mm (±0,2mm)



Magnification factor:**
30X

Optical Resolution:
3,5 µm

Field of View (FOV):
2mm

**These parameters may vary when the lenses are used with different smartphone models.*

*** Magnification obtained on a 6.4" screen smartphone.*



DIPLE®

is a SMO's
Registered Trademark

Made in Italy by
SmartMicroOptics srl

info@smartmicrooptics.com
www.smartmicrooptics.com



WARNING:


CHOKING HAZARD - Small parts.

Not for children under 3 years



DIPLE Lux:

The Diple Lux objective lens is connected to the handle. Inside the handle there is the circuit for the lighting of the sample



Batteries:




Two CR2032 (1.5V) are inside.

For changing/mounting the batteries use a screwdriver to release the three screws on the plastic shell.

Technical Specifications*:

Working distance (WD): 2.8mm (±0,2mm)



Magnification factor:**
30X

Optical Resolution:
3,5 µm

Field of View (FOV):
2mm

**These parameters may vary when the lenses are used with different smartphone models.*

*** Magnification obtained on a 6.4" screen smartphone.*



DIPLE®

is a SMO's
Registered Trademark

Made in Italy by
SmartMicroOptics srl

info@smartmicrooptics.com
www.smartmicrooptics.com



WARNING:

CHOKING HAZARD - Small parts.

Not for children under 3 years